

XML Lab Considerations

1. The database name is ISERIES
2. Page 34 to call the program **qdbxm.qzxadm** use the Qsh interface as indicated on the Note box on the same page. To do this use a 5250 emulation session defined on the desktop as As400ws.
2. Page 53 to call the program **dxsamples.shdx** use a 5250 session to do it. And the name of the database is not I400WS it is ISERIES, so replace the parameter I400WS for ISERIES.
3. Page 62 to call the program **dxsamples.chgccsid** use a 5250 session to do it. Use instead the program **dxsamples.cf** and note that the 5th parameter is blank. (")
4. Page 63 do not use Internet Explorer to visualize the file. Use Wordpad instead.
5. Page 72 to call the program **dxsamples.genx** use a 5250 session to do it.
6. On page 72 step 6 use wordpad instead.
7. On page 73 to the program call **dxsamples.genx** use a 5250 session to do it.
8. On page 73 step 4 use wordpad instead.
9. On page 77, to execute the programs **dxsamples.genx** and **dxsamples.chgccsid** use the 5250 session to do it.
10. On page 79, to execute the programs **dxsamples.genx** and **dxsamples.chgccsid** use the 5250 session to do it.